=> DOMINION
WAINARY PROTEINS
TEST FOR PROTEINS IN WAINE
and a state there methods as Identifiance
probeins in wrine (in the Laboratory).
property in wither the tobolivery.
ANSWER
(1) Heat Coagulation Test.
(2) Sulpho salityclic and TESE (SSA)
(3) Urine dip Strip Test.
Question 2 -> List three Causes of probeins in
wine (proteinuria)/bisease Conditions that Guess
proteinuria.
ANSWER
(1) Urinary Track injection
2 Multiple myeoloma
3 kidney Stone.

EXPERITION EXPENSES SUIPHOR ALM TO APPARATIONS: BEAR PLANTS SUIP PRINCIPLES Charles Turbidity Clic of Turbidity Capacity	salicyclic acid test test for protein is test for protein is to best testube, testube is to besalicyclic acid Co ck manual alke Int of Sample acid (SSA) and Che sy denotes the present ally albumin.	In wine. Tack posteur protein  BA and B and  Bands of sulphosticy  ack for protein
TEST	OB SEAUATION	
Iml of Sample At 3mls of SSA	Precipitation occurs	probeins are
	the protection does	Absence of
3mb of SSA	I not occur	Proteins.
DV8 PAST QUESTIONS		
(1) State/List three methods of identifying protein in Unite Answer > Already answered above.		
Anguar -	ruckey annuared	
Answer >	allos Do makerin	110111 000
Answer >	Causes of protessura at lead to protessura above.	110111 000

- Page 14 PAECIPITATION BY HEAT COAGULATION E : HEAT COAGULATION TEST AIM : To test for protein in wine APPARATUS: Test tube, test tube rack, water bath, pasteur Pipelle. BEAGENT: Accesiocoldacid, Sample A and B PRINCIPLE: Etack manual PROCEDURE: Hate 5ml of Sample A and Sample B in a best tube fornd heat for 5 minutes on a water bath and Coloquiation would occur if there is protein. Acidiff both samples A and B with 3-5 drops to actic acid, adding the acid drop by drop to the hot Solution. Phosphales pre capitate will disapper under these Conditions depending on the present or absent of protein in the Solution. BESULT OBSERVATION MEGRENCE Uml of Sample A Coagulation does not problem is + 5 mms heat oder absent precipitable does not protein a abent + 3 dies of actic ood in drops disappear ( lm l op sample & protein is present Coagulation occurs + 5 drops of protein is present precipitate disapeas


PAST QUESTIONS (1) Mention one Characteristic of Bence-Jopes protein and name the disease it is a bomarker ANSWER The Characteristic preperty of Beng-Jones protein is that & Coogulates between 45-60°C and reassolus above 60°C. NOTE > Normal proteins Coagulates at 68°C & don't redspolve It is a biomarker for Multiple Myeloma (2) How is Multiple mycoloma diagnosed in the lab. ANSWER The presence of Bence yours preten in work is an Important factor in the diagnosis of multiple mycelling in the Laboratory Bence Jones pretern Coagulates between 45-60°C and refussolus above 60°C. Take up the Sample to be tropped (a little 9 wantsty 5ml) and heat for about 5 minutes, if the Sample Coagulates between 45-600 and redissolves above age than there is multiple mycolama

Pag 10 EXP 2 > PURIFICATION OF PROTEINS AIM To purify or separate protein in a given sample APPARATUS: Spatula, Test tuke, test tuke rack beaker, pasker Pipette. Salt used for the Experiment.

AMMORIUM Sulphabe, Sample A and B PRINCIPLE: Check manual. Also Stated below. PROCEDURE: To I'm of test solution, add little by little ammonium sulphake and observe; Continue adding until the Solution precipitate. PAST QUESTION (1) State Clearly Salting- in and Salting Out (Principle ANSINER when a Crystal of a Soluble Salt (6.9 Nach) are added little by little to a protein Solution, the first few Crystals of the Salt go into the Solution, being dissolved by water of hydration attached to the protein Molecule. This phenomenon & Called Salting In On further addition of the Selts the water of hydron is gradually lost to increasing Salt mokeule. Consequently, the protein begins to precipitate. At Saturation point, all the protein will precipitate.
This is Could Salting Out

EXPERIMENT 3 LE: Gastric acidity determination ALM: To determine the gostic acidity Content of the body. Pasteur APPARATUS: Conical flask, pipette, beaker chock manual for REAGENT: Thymol blue indicator pronciple. PROCEDURE: Take Ambs of Sample A and Sample B in your Conical flask. TO Sample A - Add 3 drops of Khymol blue indicator - Titrate with 0.1 Al of Nach. One you get the first Colour Change to yellow note the amount of Conc baot - That becomes your free acidity, Continue until you get permanent blue Colour, amount again - The first titre value and the second bitre value Summed up is the total acidity, - Do Same for Sample B. pht is defined as the measure of the degree of a cidity or alkalishy in a Solution.

@ what is pH? In Chemistry, pH historically denoting Gotential of hydrogen" (or "Grower of hydrogen") is a scale used to specify the acidity or basicity of an aqueous solution 3) What is the name of the indicator of gastric acid analysis? ANSWER Thymol blue. NOTE SITE has a different ph ranges that is why it is wed (4) State two acids being determined in gastric Juie analysis ANSWER 1) free acid (1) Total acid (3) what is the signifiance of gostric juice analysis ANSWER It is used to determine pelptic ulcer disease OTHER NOTABLE POINTS 6 Constituents of gastric juillacid Free HCL 2 weak acid such as phosphate, protein and organicais 3) Organic acid such as acetic acid, lactic acid.




page 11 EXPENIMENT 4 Nin hydrin Test To determine or identify a-amino acid APPARATUS: Beaker, test tuke rack, pasteur prette
REAGENT: Ninhydria reagent. Check manual
principle REAGENT: PROCEDURE: TO Hall of Sample A, add 2 drops of Non hydrin reagent. Then heat watch for the purple Colouration. Do same for Sample B, C and & respectively. Record your result. RESULT INFERENCE OBSERVATION Int of Sample A+ Purple Colouration present of amino 2 drops of Norhydin of and acid acid reagent + 2 mins heat presence of amino I'ml of Sample Bt purple Colouration acid 2 drops of Ninhydrin reagent + 2 mins hook Int of Semple C+ Large of Semple C+ Nonhydrin reagent + Yellow Colouration 2 mins heat

DIS CUSSIEN Non hydran test is a unwested method that is wed to assay or determine amino acid, frimary amines, Secondary aminus respectfully of-amino acide. BIGHTRAHCAL IMPORTANCE It is used for finger print detection. => warme two among acred that does not give purple Colour in unpydrin Lest. ANSWER 1 Profine 8 hydroxyprotin which gives a yellow Colour. @ Asparagine which gross a brown Colour. => Another name Nonhydrin reagent ANSWEA triketohydrindene hydrate: => The appearance of purple Colour on Northydron tox Indicates? ANSWER The present of amono acids > State Ninhydrin principle.

Page 11 TITLE: BOUTCE TEST APPARATUS: Test tube, posteur pipette, Test tube rock, Beaker. AFAGFOT: Buiret reagent, Normal Salme, plasma. PROCEDURE: To I'm of plasma, add that of Burret reagent - Allow to stand at room temperature for 15 minutes, he cord your observation.
- Carryout the first on 1/2, 1/4, 1/8, dilution Of your test solution, wing Bourne Drum albalin dilution. - Allow Colour to develop for 5 minutes and record Colour intensities at 550mm wave length in a Calorimoter. PAINCIPLE: Compounds Confaming two or more perside bonds give a Characters the purple Colour when tracked thith ditute Copper Sulphate in alkalise Solution The name of the test Comes from the Compound burret which gives a typically positive reaction.
The colour is apparently due to the co-ordination Complex of the Copper atom and 4 nitrogen abouts, two from each of the two peptide chairs.

NOTE -> Alkaline Solution is used morder not to benobe proteins

DISCUSSION
what is protein?
protein are polymers of amino acids joined together
by perfide bond.
what is are amno acids?
Amino acid care organic Compounds that has bett
NH2 and Coo group.
What Colour is gotten from Buret Test?
Violet Colour
RESULT
Serum Neat /2 /4 1/8
Plasma (m) m)
Salvie - Im/ -> Im/ -> Im/
Bivet 4ml 4ml 4ml 4ml
Reading   ++++   ++   +
@ Lobat box of the strong
Derum. Ger you use to identify protein in
Birret Test Answer
OTUCE TOOK
1 Lost two anti-Coaquiants used in the Laboratory
ANSWER
1) Marchael (4) CM
EDTA Contrate (They all remove)



3 DIFF FRENCE BETWEEN PLASMA AND SERVIN SERUM (1) This is kno liquid that this is the liquid that remains when anticoagulant remouns after bho abbing is added to prevent clothing of blood less volume Computed to (2) Consolls of 55% of the total blood volume. dosma Lacks Polarinogen, (3) Conteurs Fibringen (4) It consists of serum and It is the part of blood without clothing Packer Clothing factor. 9 It needs Entitogulant Anticogulant of not before it can be abbained pood before it can be Ototewned from blood Sample (6) List the plasma proteins on the basis of their solvebility in Salt Solutions-ANSW ER Fibrinogen (11) Globulin (III) Albumin (7) what are the Conditions that (nereases plasma Moterns in the body,

ANSWEA Dehydra Loon Doarrhoea Multiple myedoma Posporatory disorder 5 vomitting exastively Respiratory disorder Hemoly869 Inflamatory disorder. Conditions that becreases Plasma Proteins Malnutrition & Cas seen in kwashiotor) Cirhossis of the Liver Malabsorption Sovere burns Nophritic Syndrama.